



GSVMX-S



GSVMX-F

Features

- Variable speed of 0 to 3,000rpm for gentle shaking to rapid mixing of samples
- Touch or continuous operation mode
- Includes a standard top head for mixing a single test tube or other small vessel up to 30mm in diameter
- Electronic speed control maintains constant speed even during changes in load
- Steel base with specially designed vacuum suction feet provides stability and dampens vibration
- Robust aluminum-cast construction

Specification

Model	GSVMX-S	GSVMX-F
Speed range	0-3000rpm	3000rpm
Speed display	Scale	/
Run type	Touch operation and Continuous	
Mixing motion	Orbital	
Orbital diameter	4mm	
Motor type	Shaded pole motor	
Motor rating input	58W	
Motor rating output	10W	
Permissible ambient temperature	5-40°C	
Permissible relative humidity	80%RH	
Protection class	IP21	
Power Supply	110-120V/220-230V,50/60Hz	
Consumption	60W	
External Size (W×D×H)(mm)	130×127×160	
Shipping Size (W×D×H)(mm)	300×250×230	
N.W./G.W. (kg)	2/3.5	