

HPLC COLUMN STORAGE CABINET

 **PRECISE** |  **SAFE** |  **SECURE**



*Terms and Conditions Apply

SPECIFICATIONS

Product Name	HPLC Column Storage Cabinet
Construction	Solid CRC Steel duly powder coated with drawer handles
Column Length Capacity	Upto 30 cm
Columns Capacity	50 HPLC Columns
Drawer Type	Multiple Drawers with Column Holding Racks
Locking System	Key Lockable for Extra Safety
Finish	Powder Coated
Application	Storage and Protection of HPLC Columns
Color	Light Grey
Usage	Laboratories, Research Institutes, QC/QA Departments